

# Clean-in-Place Analysis Applications [OMA-300-CIP]

Applied Analytics Application Note No. DS-001F — Revised 24 June 2013



## Application Summary

Analyte: Trace impurities in a cleaning solvent (e.g. water, methanol)

Detector: OMA-300 Process Analyzer (using UV-Vis / SW-NIR dispersive absorbance spectrophotometer)

Process Stream: Cleaning solvent

## Introduction

Clean-in-Place (CIP) is a method commonly used to clean the interior surfaces of production equipment. In the process, cleaning solvents such as methanol or water are circulated to sanitize walls, pipes, tanks, etc. To satisfy regulations, the wash must effectively remove all active ingredients from previous production run.

CIP is a stage of significant downtime; production will be stalled until the process is complete. Plants that still use a set duration or depend solely on swab tests are suffering unnecessary downtime and costs. Allowing CIP to run longer than it needs to is an oversight that wastes production time and CIP solvents.

Using an OMA system, the wash cycle is validated in real time. An active ingredient (the substance which is being eradicated) will typically show distinct absorbance curves in the UV-Vis or SW-NIR ranges; the OMA monitors this absorbance to measure the current concentration of active ingredient X in the effluent cleaning solvent (rinse water). When the threshold concentration is reached in the continuously drawn rinse water sample, the OMA immediately notifies the user or external system. This functionality reduces unneeded CIP downtime and removes dependency on labor-intensive swab tests.

## OMA Benefits

- » Automates real-time validation of Clean-in-Place wash process
- » Faster response than alternative technologies (e.g. ion mobility)
- » Rugged, solid state build with no moving parts — modern design for low maintenance
- » Long-lifespan xenon light source (avg. 5 years)
- » Specialized sanitary flange sample flow cell
- » Optional remote access

# Clean-in-Place Analysis Applications [OMA-300-CIP]

Applied Analytics Application Note No. DS-001F — Revised 24 June 2013


## Sanitary Flange Flow Cell

The OMA for Clean-in-Place analysis uses a special low cell design with sanitary quick-clamp flanges. The flow cell is constructed from 316L stainless steel with a path length dependent upon the measurement application.



## Absorbance Curves of Example Active Ingredients

These spectra visualize how the OMA correlates the concentration of a unique active ingredient to the level of its absorbance curve in the sample.




# Clean-in-Place Analysis Applications [OMA-300-CIP]

Applied Analytics Application Note No. DS-001F — Revised 24 June 2013


## Clean-in-Place Application Schematic

Some of the recycled solvent is diverted to a small bypass flowing through the sanitary flange flow cell. After reaching a steady state concentration in the flow cell, a new batch of cleaning solvent is loaded and the process is repeated until the impurity concentration reaches the next preset threshold value.



## Clean-in-Place Timing Diagram

The concentration of the active ingredient in the solvent increases until it reaches a plateau, signaling the end of the wash. The 1st derivative of the concentration with time ( $dc/dt$ ) approaches zero when the process reaches steady state.



# Clean-in-Place Analysis Applications [OMA-300-CIP]

Applied Analytics Application Note No. DS-001F — Revised 24 June 2013

## Further Reading

| Subject  | Location  |
|--|---|
| OMA-300 Process Analyzer<br>Data sheet         | <a href="http://www.a-a-inc.com/documents/AA_DS001A_OMA300.pdf">http://www.a-a-inc.com/documents/AA_DS001A_OMA300.pdf</a>                 |
| Advantage of Collateral Data<br>Technical Note | <a href="http://www.a-a-inc.com/documents/AA_TN-202_CollateralData.pdf">http://www.a-a-inc.com/documents/AA_TN-202_CollateralData.pdf</a> |



is a registered trademark of Applied Analytics Group BV. | [www.a-a-inc.com](http://www.a-a-inc.com)

### Headquarters + Manufacturing

Applied Analytics, Inc.  
Burlington, MA | [sales@a-a-inc.com](mailto:sales@a-a-inc.com)

### North America Sales

Applied Analytics North America, Ltd.  
Houston, TX | [sales@appliedanalytics.us](mailto:sales@appliedanalytics.us)

### Europe Sales

Applied Analytics Europe, SpA  
Milan, Italy | [sales@appliedanalytics.eu](mailto:sales@appliedanalytics.eu)

### Asia Pacific Sales

Applied Analytics Asia Pte. Ltd.  
Singapore | [sales@appliedanalytics.com.sg](mailto:sales@appliedanalytics.com.sg)

### Middle East Sales

Applied Analytics Middle East (FZE)  
Sharjah, UAE | [sales@appliedanalytics.ae](mailto:sales@appliedanalytics.ae)

### Brazil Sales

Applied Analytics do Brasil  
Rio de Janeiro, Brazil | [sales@aadbl.com.br](mailto:sales@aadbl.com.br)

### India Sales

Applied Analytics (India) Pte. Ltd.  
Mumbai, India | [sales@appliedanalytics.in](mailto:sales@appliedanalytics.in)

© 2013 Applied Analytics Group BV. Products or references stated may be trademarks or registered trademarks of their respective owners. All rights reserved. We reserve the right to make technical changes or modify this document without prior notice. Regarding purchase orders, agreed-upon details shall prevail.